

# Logon from CME Globex to Client System

This topic contains the message specification for the Logon (tag 35-MsgType=A) message sent by CME Globex in response to Logon (tag 35-MsgType=A) message received from client system.

Tag	Name	FIX Type	Enumeration	Req	Description
<a href="#">Standard Header - CME Globex to Client System</a>					
35	MsgType	String (2)	A=Logon	Y	Header tag identifying message type.
108	HeartBtInt	Int(3)	30 - 60 inclusive	Y	Heartbeat interval in seconds. CME Globex will echo back the heartbeat interval submitted in the customer logon message.
141	ResetSeqNumFlag	Boolean(1)	Y=Yes, reset sequence numbers N=No	C	CME Globex echoes back the ResetSeqNumFlag indicator from the customer logon message.
1603	ApplicationSystemName	String (30)	Echoes the value sent from the client system.	Y*	Identifies system generating the message. <b>Note:</b> this tag is returned if in the header for the <a href="#">Logon from CME Globex to Client System</a> message only.
1604	TradingSystemVersion	String (10)	Echoes the value sent from the client system.	Y*	Identifies system generating the message. <b>Note:</b> this tag is returned if in the header for the <a href="#">Logon from CME Globex to Client System</a> message only.
1605	ApplicationSystemVendor	String (10)	Echoes the value sent from the client system.	Y*	Identifies system generating the message. <b>Note:</b> this tag is returned if in the header for the <a href="#">Logon from CME Globex to Client System</a> message only.
<a href="#">Standard Message Trailer</a>					
					End of message.

Y: Required by FIX protocol, Y\*: Required by CME Globex (not by FIX protocol), N: Not Required, C: Conditionally